

Practice baby for nursing training

Item no.: BA186 Nora

Factsheet



*Experts in medical education
since 1950*



ERLER ZIMMER

Erlers-Zimmer GmbH & Co. KG · Hauptstr. 27 · D-77886 Lauf
Tel.: (+49) 7841/67191-0 · info@erler-zimmer.de · www.erler-zimmer.de

Practice baby for nursing training

Nurses face very special challenges in infant nursing. This is why careful training in this area is so extremely important. With the Nora nursing baby, you can prepare your trainees ideally for this demanding task.

Our new nursing manikins offer practice opportunities for all important nursing measures on the infant.

The following procedures are possible on the manikin:

- Gastric tube (placement, feeding, irrigation and suctioning)
- Tracheotomy (insertion and suctioning)
- Gastrostomy / PEG (care)
- Suprapubic catheter (changing, care)
- Bladder catheterization (insertion, placement and care)
- Ileostomy, colostomy, double-barreled stoma (basic care only)
- Inserting rectal thermometers
- Basic care measures such as washing, wrapping, etc.

