

# Erwachsenen-Ganzkörperpuppe

für Röntgen, CT, US ohne Muskeln

## Adult Full Body manikin

for X-Ray, CT, US without muscles



Merkblatt / Factsheet

Art. Nr. / Item No. FLB03

## Erwachsenen-Ganzkörperpuppe

Dieses Ganzkörper-Phantom für Erwachsene ist ein Röntgen/CT- und US-kompatibles Trainingsprodukt. Die Hauptanwendung dieses Phantoms ist die praktische Erfahrung mit diagnostischen Bildgebungsverfahren und die visuelle Bewertung, um optimale Scanbedingungen zu finden.

Das Modell dieses Phantoms basiert auf einem erwachsenen männlichen Körper, der etwa 70 kg wiegt und 172 cm groß ist.

Die Knochen haben eine realistische dreischichtige Struktur mit innerer Porosität, die je nach den Anforderungen des jeweiligen Projekts angepasst werden kann. Auf Wunsch mit verschiedenen Pathologien (Läsion, Tumor, Infektion, Anomalie) anpassbar. Dieses Positionierungsphantom kann in einer durchsichtigen oder hautfarbenen Version bestellt werden.

### Anatomie:

- Vollständiger menschlicher Körper (10 Teile)
- Realistischer Körper und muskelgewebeähnliches Material
- Drehbare Schultern (360° rundum und 180° seitlich)
- Bewegliche Hüft-, Knie- und Ellbogengelenke
- Kopf, Torso und Gliedmaßen abnehmbar
- Erwachsener Kopf
- Erwachsener Torso
- Arme des Erwachsenen
- Erwachsene Beine

## Adult Full Body manikin

This Adult Full Body Phantom is an X-Ray/CT and US-compatible training product. The primary application of this phantom is hands-on experience with diagnostic imaging techniques and visual evaluation to find optimal scanning conditions.

The model of this phantom is based on an adult male body that weighs about 70 kg (154 lb) and measures 172 cm (69 in) in height.

The bones have a realistic three-layered structure with inner porosity, which can be adjusted according to the requirement of the particular project. Customizable with different pathologies (lesion, tumor, infection, abnormality) upon request. This positioning phantom can be ordered in a translucent or skin-colored version.

### Anatomy:

- Full Human Body (10 Parts)
- Realistic Body and Muscle Tissue-Mimicking Material
- Rotatable Shoulders (360° Around and 180° Sideways)
- Movable Hip, Knee, and Elbow Joints
- Detachable Head, Torso, and Limbs
- Adult Head
- Adult Torso
- Adult Arms
- Adult Legs

