## **Crico Chris - Cricothyrotomy**

**Factsheet** 

**Task Trainer** 





## Crico Chris - Cricothyrotomy Task Trainer

Cricothyrotomy places high demands on the skills of the practitioner and can have serious consequences for the patient if not practiced properly. Crico Chris enables regular training of lifelike cricothyrotomies without risk. The handy trainer has a larynx with interchangeable conic ligament and skin. Both the skin and Llg. Conicum are designed so that both can be gradually rotated 360° and allow a variety of exercises until replacement is required. A special feature of Crico Chris is the blood-filled pouch placed between the skin and ligament, which allows blood to escape during the incision and makes further work more difficult. In addition to lifelike anatomy and a realistic feel, Crico Chris offers an excellent price/performance ratio and high cost-effectiveness.

## Scope of delivery:

Basic body

- Laryngeal cartilage
- 2x Larynx ligament layer
- Throat skin
- Airway
- 3x Artificial blood bag

## Spare parts / consumables:

Item no.: 7660A Skin

**Item no.: 7660B** Ligaments layer, set of 3 **Item no.: 7660C** Artificial blood bag, set of 3

**Item no.: 7660D** Laryngeal cartilage

Item no.: 7660E Airway







