

MILITARY WORKING DOG SIMULATOR K9-MARSHAL

Item no.: VET4700








The “**Marshal**” military working dog simulator represents a **Belgian Malinois** and corresponds in size and weight (30 kg) to a real dog and has **movable joints**.

Military working dogs are used in many areas by the military and police and, like their human comrades, can be injured in the line of duty. **Marshal helps train vets, K9 Combat First Responders and dog handlers.**

With Marshal you can present a variety of realistic scenarios during training.

Made of silicone with internal skeleton. Body cover made of robust material so that Marshal can be exposed to adverse environmental influences such as moisture and dirt without being damaged.





Marshal offers the following features:

- 1 Lacerated paw for bandage training
- 2 Manual femoral pulse control
- 3 Bullet wound for wound packing – Stop the bleed
- 4 Basic cardiopulmonary resuscitation (CPR) techniques (without chest lift)
- 5 2-layer exchangeable fur (with robust zippers and sewing)
- 6 Moveable head
- 7 Movable jaw for muzzle training / dog bite training
- 8 Open nostrils
- 9 Endotracheal intubation (without chest lift)
- 10 Flexible tongue
- 11 Compartment with Bluetooth speaker (different barking, whining and crying sounds available)
- 12 Movable joints
- 13 Exploded bleeding stump for tourniquet training (non-exchangeable)

Dog vest and goggles not included.



EXPERTS IN MEDICAL EDUCATION

Erlor-Zimmer GmbH & Co.KG

Hauptstraße 27 · 77886 Lauf · Germany

T (+49) 7841 / 67191-0 · F (+49) 7841 / 67191-99

info@erler-zimmer.de

www.erler-zimmer.de

Follow us!

